

## Paranasal Sinus Model



### Characteristics

- illustration of the bone structures of the paranasal sinus system (for example: ethmoidal cells, maxillary caves, frontal sinus, sphenoidal sinus)
- reproduction of soft tissue (for example: skin, dura, nervus opticus, aorta carotis interna, aorta ethmoidalis anterior/ posterior)
- representation of special disease pictures and coloured accentuations are possible
- the model consists of a mineralogical powder material

### Fields of application

- education for medical specialists
- courses of instruction / training
- visualisation

